

# TT4 Rapid Quantitative Test

Catalog No.BT2210

## INTENDED USE

The Biotime TT4 Rapid Quantitative Test is intended to quantify the concentration of total thyroxine in human serum or plasma on Biotime FIA Analyzer by fluorescent immunoassay. The test is used as an aid detection of thyroid gland disease.

-Fluorescence immunoassay.

-Thyroid gland disease test.

For in vitro diagnostic use only. For professional use only.

## INTRODUCTION

T4 is produced by the thyroid gland. Adults bodies normally produce about 80-100µg per day. 99.96% of the T4 in the blood exists as protein-binding form, 90% of which binds to thyroid hormone binding globulin (TBG).

Quantification of TT4 helps to evaluate thyroid functions. Blood TT4 levels elevate in patients with hyperthyroidism and decrease in patients with hypothyroidism. In patients with declined thyroid functions, TT4 levels are also lower than normal conditions.

TT4 measurements is related to the concentrations of its carrier proteins. Thus evaluation of T4 binding capacity to its carrier proteins is necessary. TT4 diagnostic results should also be combined with other laboratory test or clinical examines to evaluate thyroid functions<sup>[1-3]</sup>.

## PRINCIPLE

This reagent is based on fluorescent lateral flow immunoassay competition method.

While the sample and the buffer are mixed, T4 in specimen and mouse anti-T4 monoclonal antibody labeled with fluorescent microsphere(contain Europium) form a reaction complex. During lateral flow, free labeled antibodies move along with the nitrocellulose membrane to a detection line (T-line:coated with T4 antigen). Free labeled antibodies are captured and gives fluorescent signal upon stimulation. Thus, the fluorescent signals are negatively correlated with the concentrations of TT4 in human plasma or serum.

The fluorescent signal will be quantified and calculated according to the calibration curve (provided with the reagents) to represent concentration of TT4 in specimen.

## PRECAUTIONS

1. This reagent is used for in vitro diagnostic use only, please do not use expired products.
2. All blood samples (including the remaining samples after testing), used reagents and waste should be treated as infectious materials.
3. The reagent is for one-time use. Once the pouch is opened, it should be used within 30 minutes to avoid failure caused by the moisture absorption.
4. While using the test cartridge and instruments, vibration and electromagnetic environment should be avoided.
5. Lot number of buffers and test cartridges must be matched.
6. Do not insert the cartridges that are contaminated with blood or other liquids on the surface. It may cause damages to the instrument.

## MATERIAL

### Material Provided

1. Test cartridge 25 tests/kit
2. Detection buffer 25 tubes/kit
3. SD card 1 piece/kit
4. Instructions for use 1 copy/kit

### Material Required But Not Provided

1. Biotime FIA Analyzer
2. Transfer Pipette Set (5~50uL and 10~100uL size)
3. Specimen collection containers
4. Timer

## STORAGE AND STABILITY

1. Store the detection buffer at 2-8°C, the shelf life is 24 months.
2. Store the test cartridge at 2-30°C, the shelf life is 24 months.
3. Test Cartridge should be used within 30 minutes after opening the pouch.

## SPECIMEN COLLECTION AND PREPARATION

1. The test can be performed with plasma or serum.
2. The collection of sample container should be heparin anti-coagulant tube for plasma or pro-coagulant tube for serum.
3. The collection of the sample: the venipuncture for blood collection method referring to the National Clinical Laboratory Procedures, if the sample can't be detected timely, it must be stored in refrigerator at 2-8°C for no more than 7 days, or at -20°C for no more than 6 months. Samples should be recovered to the room temperature before tests.
4. Separate the plasma or serum from blood as soon as possible to avoid hemolysis.

## TEST PROCEDURE

Please refer to the operation manual of Biotime FIA Analyzers for details.

The test should be operated at room temperature (~25°C) .

### Step 1: Preparation

Check/insert SD card into the instrument.

Take out one tube of buffer from refrigerator and balance it to room temperature.

### Step 2: Sampling

Take 40µL of plasma or serum with a transfer pipette and add it to the buffer tube.

### Step 3: Mixing

Mix well the specimen with buffer by tapping or inverting the tube.

### Step 4: Loading

Take 80µL of sample mixture and load it into the well of the test cartridge.

Note: Step 2 to step 4 should be completed within 1 minute to ensure the accuracy of the test results.

### Step 5: Testing

Ensure that there are no air bubbles. Immediately insert the test cartridge into analyzer and incubate for 15 minutes.

NOTE: Please refer to the operation manual of a specific model of the analyzer for details.

## REFERENCE INTERVAL

Normal Reference Value:58.0-140.0nmol/L

Note: Individual reference range is suggested to be established for each laboratory..

Unit conversion formula: ng/mL\*1.287=nmol/L

## LIMITATIONS OF PROCEDURE

1. The test sample should be plasma or serum.
2. Human anti-mouse antibody (HAMA) may be present in patients who have received immunotherapy with a murine monoclonal antibody. This kit has been specially designed to minimize the effect of these antibodies on the test results. However, the test result must be carefully evaluated when patients are known to have these antibodies<sup>[4-5]</sup>.
3. Other factors also can induce the false results, include the technology, operational error and other sample factors.

## PERFORMANCE CHARACTERISTICS

### Accuracy

TT4 control materials with two different concentrations were tested by every lot of Test Cartridges, and the deviations are within ±15.0%.

### Assay Range and Detection Limit

**Assay Range:** 10.0-320.0nmol/L

**The Lowest Detection Limit:** 10.0nmol/L

### Linearity

A serial concentration of TT4 reference materials at 20.0-320.0nmol/L were tested,and the correlation coefficient (R) is ≥0.9900.

### Precision

#### Intra-Lot Precision

Intra-lot precision was determined by testing of TT4 reference materials using 10 test cartridges from the same lot. The C.V.is ≤ 15.0%.

#### Inter-Lot Precision

Inter-lot precision was determined by testing of TT4 reference materials using 30 test cartridges from 3 consecutive batches randomly (10 test cartridges from each lot). The C.V.is ≤ 20.0%.

### Specificity:

The concentration of TT4 is not greater than 15nmol/L when the concentration of TT3 is 500.0ng/mL.

The concentration of TT4 is not greater than 15nmol/L when the concentration of rT3 is 50.0ng/mL.

## SYMBOLS

Symbol	Description	Symbol	Description
	Catalogue number		In vitro diagnostic medical device
	Batch code		Consult instructions for use
	Date of manufacture		Keep dry
	Use-by date		Keep away from sunlight
	Manufacturer		Store at 2-8°C


	Do not re-use		Store at 2-30°C
	Authorized representative in the European		CE mark

**BIBLIOGRAPHY OF SUGGESTED READING**

1. Li Guixiang.Detection of TT3, TT4, TSH, FT3 and FT4 and judgment of thyroid function[J]. Chinese Pharmaceutical Journal, 2003(9):81.
2. Liappis N,Schlebusch H,Mallmann R.Comparison of thyrotropin (TSH),thyroxine (TT4)and triiodothyronine (TT3) concentration in capillary and venous blood serum in children and adolescents.Monatsschr Kinderheilkd[J].1991 Nov;139(11):765-7.
3. Liappis N,Schlebusch H,Knapp M.Reference values for the concentration of free triiodothyronine (FT3), triiodothyronine (TT3), free thyroxine (FT4), thyroxine (TT4), thyrotropin (TSH) and thyroxine-binding globulin (TBG) in umbilical cord blood. Method: luminescence-enhanced enzyme immunoassay.Klin Padiatr.1992 Jan-Feb;204(1):34-6.
4. Hansen JH, et al. HAMA Interference with Murine Monoclonal Antibody-Based Immunoassays[J].J of Clin Immunoassay,1993,16:294-299.
5. Levinson SS.The Nature of Heterophilic Antibodies and the Role in Immunoassay Interference[J].J of Clin Immunoassay,1992,15:108-114.



Xiamen Biotime Biotechnology Co., Ltd.  
 Address: 2F/3F/4F/5F, No.188, Pingcheng South Road, Haicang Street,  
 Haicang District, Xiamen City, Fujian Province, 361026, P. R. China.



Lotus NL B. V.  
 Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherlands.  
 Tel: +31644168999

Version: A/09  
 Issuing date: 2022-07-01